

Antibodies Product Datasheet

Product Name:	ARRDC2 Antibody
Catalog Number:	C12060
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARRDC2. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	ARRDC2
Protein Name:	Arrestin domain-containing protein 2
Other Name:	ARRDC2; PP2703; Arrestin domain-containing protein 2
Human Swiss Prot No.:	Q8TBH0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARRDC2 Polyclonal Antibody detects endogenous levels of ARRDC2 protein.
Storage and Stability:	-20°C/1 year