

Antibodies Product Datasheet

Product Name:	AS250 Antibody
Catalog Number:	C12064
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AS250. AA range:641-690
Concentration:	1 mg/ml
Gene Name:	RALGAPA2
Protein Name:	Ral GTPase-activating protein subunit alpha-2
Other Name:	RALGAPA2; C20orf74; KIAA1272; Ral GTPase-activating protein subunit alpha-2; 250 kDa substrate of Akt; AS250; p220
Human Swiss Prot No.:	Q2PPJ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AS250 Polyclonal Antibody detects endogenous levels of AS250 protein.
Storage and Stability:	-20°C/1 year