

Antibodies Product Datasheet

Product Name:	CYB561D1 Antibody
Catalog Number:	C12081
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CYB561D1. AA range:180-229
Concentration:	1 mg/ml
Gene Name:	CYB561D1
Protein Name:	Cytochrome b561 domain-containing protein 1
Other Name:	CYB561D1; Cytochrome b561 domain-containing protein 1
Human Swiss Prot No.:	Q8N8Q1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYB561D1 Polyclonal Antibody detects endogenous levels of CYB561D1 protein.
Storage and Stability:	-20°C/1 year