

Antibodies Product Datasheet

Product Name:	CYB561D2 Antibody
Catalog Number:	C12082
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C56D2. AA range:55-104
Concentration:	1 mg/ml
Gene Name:	CYB561D2
Protein Name:	Cytochrome b561 domain-containing protein 2
Other Name:	CYB561D2; 101F6; LUCA12.2; Cytochrome b561 domain-containing protein 2; Putative tumor suppressor protein 101F6
Human Swiss Prot No.:	O14569
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYB561D2 Polyclonal Antibody detects endogenous levels of CYB561D2 protein.
Storage and Stability:	-20°C/1 year