

Antibodies Product Datasheet

Product Name:	Cables1 Antibody
Catalog Number:	C12084
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Ik3-1. AA range:561-610
Concentration:	1 mg/ml
Gene Name:	CABLES1
Protein Name:	CDK5 and ABL1 enzyme substrate 1
Other Name:	CABLES1; CABLES; CDK5 and ABL1 enzyme substrate 1; Interactor with CDK3 1; Ik3-1
Human Swiss Prot No.:	Q8TDN4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cables1 Polyclonal Antibody detects endogenous levels of Cables1 protein.
Storage and Stability:	-20°C/1 year