

Antibodies Product Datasheet

Product Name:	K-cadherin Antibody
Catalog Number:	C12086
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDH10. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	CDH10
Protein Name:	Cadherin-10
Other Name:	CDH10; Cadherin-10; T2-cadherin
Human Swiss Prot No.:	Q9Y6N8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	K-cadherin Polyclonal Antibody detects endogenous levels of K-cadherin protein.
Storage and Stability:	-20°C/1 year