

Antibodies Product Datasheet

Product Name:	Cadherin-23 Antibody
Catalog Number:	C12097
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDH23. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	CDH23
Protein Name:	Cadherin-23
Other Name:	CDH23; KIAA1774; KIAA1812; Cadherin-23; Otocadherin
Human Swiss Prot No.:	Q9H251
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cadherin-23 Polyclonal Antibody detects endogenous levels of Cadherin-23 protein.
Storage and Stability:	-20°C/1 year