

Antibodies Product Datasheet

Product Name:	Cadherin-24 Antibody
Catalog Number:	C12098
Reactivity:	Human, Mouse
Applications:	IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDH24. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	CDH24
Protein Name:	Cadherin-24
Other Name:	CDH24; CDH11L; Cadherin-24
Human Swiss Prot No.:	Q86UP0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cadherin-24 Polyclonal Antibody detects endogenous levels of Cadherin-24 protein.
Storage and Stability:	-20°C/1 year