

Antibodies Product Datasheet

Product Name:	IGSF4D Antibody
Catalog Number:	C12109
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CADM2. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	CADM2
Protein Name:	Cell adhesion molecule 2
Other Name:	CADM2; IGSF4D; NECL3; Cell adhesion molecule 2; Immunoglobulin superfamily member 4D; IgSF4D; Nectin-like protein 3; NECL-3; Synaptic cell adhesion molecule 2; SynCAM 2
Human Swiss Prot No.:	Q8N3J6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IGSF4D Polyclonal Antibody detects endogenous levels of IGSF4D protein.
Storage and Stability:	-20°C/1 year