

Antibodies Product Datasheet

Product Name:	IGSF4C Antibody
Catalog Number:	C12111
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CADM4. AA range:339-388
Concentration:	1 mg/ml
Gene Name:	CADM4
Protein Name:	Cell adhesion molecule 4
Other Name:	CADM4; IGSF4C; NECL4; TSLL2; Cell adhesion molecule 4; Immunoglobulin superfamily member 4C; IgSF4C; Nectin-like protein 4; NECL-4; TSLC1-like protein 2
Human Swiss Prot No.:	Q8NFZ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IGSF4C Polyclonal Antibody detects endogenous levels of IGSF4C protein.
Storage and Stability:	-20°C/1 year