

Antibodies Product Datasheet

Product Name:	Cdc50C Antibody
Catalog Number:	C12120
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TMEM30C. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	TMEM30C
Protein Name:	Cell cycle control protein 50C
Other Name:	TMEM30C; CDC50C; Cell cycle control protein 50C; Transmembrane protein 30C
Human Swiss Prot No.:	A0ZSE6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc50C Polyclonal Antibody detects endogenous levels of Cdc50C protein.
Storage and Stability:	-20°C/1 year