

Antibodies Product Datasheet

Product Name:	CD2BP2 Antibody
Catalog Number:	C12146
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD2 Tail-binding. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	CD2BP2
Protein Name:	CD2 antigen cytoplasmic tail-binding protein 2
Other Name:	CD2BP2; KIAA1178; CD2 antigen cytoplasmic tail-binding protein 2; CD2 cytoplasmic domain-binding protein 2; CD2 tail-binding protein 2
Human Swiss Prot No.:	O95400
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD2BP2 Polyclonal Antibody detects endogenous levels of CD2BP2 protein.
Storage and Stability:	-20°C/1 year