

Antibodies Product Datasheet

Product Name:	DCL-1 Antibody
Catalog Number:	C12150
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD302. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	CD302
Protein Name:	CD302 antigen
Other Name:	CD302; CLEC13A; DCL1; KIAA0022; CD302 antigen; C-type lectin BIMLEC; C-type lectin domain family 13 member A; DEC205-associated C-type lectin 1; Type I transmembrane C-type lectin receptor DCL-1; CD antigen CD302
Human Swiss Prot No.:	Q8IX05
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DCL-1 Polyclonal Antibody detects endogenous levels of DCL-1 protein.
Storage and Stability:	-20°C/1 year