

Antibodies Product Datasheet

Product Name:	Cdc27 Antibody
Catalog Number:	C12168
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human H-NUC. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	CDC27
Protein Name:	Cell division cycle protein 27 homolog
Other Name:	CDC27; ANAPC3; D0S1430E; D17S978E; Cell division cycle protein 27 homolog; Anaphase-promoting complex subunit 3; APC3; CDC27 homolog; CDC27Hs; H-NUC
Human Swiss Prot No.:	P30260
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc27 Polyclonal Antibody detects endogenous levels of Cdc27 protein.
Storage and Stability:	-20°C/1 year