

Antibodies Product Datasheet

Product Name:	CdcA4 Antibody
Catalog Number:	C12175
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDCA4. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	CDCA4
Protein Name:	Cell division cycle-associated protein 4
Other Name:	CDCA4; HEPP; Cell division cycle-associated protein 4; Hematopoietic progenitor protein
Human Swiss Prot No.:	Q9BXL8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CdcA4 Polyclonal Antibody detects endogenous levels of CdcA4 protein.
Storage and Stability:	-20°C/1 year