

Antibodies Product Datasheet

Product Name:	CDO Antibody
Catalog Number:	C12181
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDON. AA range:511-560
Concentration:	1 mg/ml
Gene Name:	CDON
Protein Name:	Cell adhesion molecule-related/down-regulated by oncogenes
Other Name:	CDON; CDO; Cell adhesion molecule-related/down-regulated by oncogenes
Human Swiss Prot No.:	Q4KMG0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CDO Polyclonal Antibody detects endogenous levels of CDO protein.
Storage and Stability:	-20°C/1 year