

Antibodies Product Datasheet

Product Name:	CDHF9 Antibody
Catalog Number:	C12182
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CELSR1. AA range:921-970
Concentration:	1 mg/ml
Gene Name:	CELSR1
Protein Name:	Cadherin EGF LAG seven-pass G-type receptor 1
Other Name:	CELSR1; CDHF9; FMI2; Cadherin EGF LAG seven-pass G-type receptor 1; Cadherin family member 9; Flamingo homolog 2; hFmi2
Human Swiss Prot No.:	Q9NYQ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CDHF9 Polyclonal Antibody detects endogenous levels of CDHF9 protein.
Storage and Stability:	-20°C/1 year