

Antibodies Product Datasheet

Product Name:	CDHF10 Antibody
Catalog Number:	C12183
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CELSR2. AA range:2781-2830
Concentration:	1 mg/ml
Gene Name:	CELSR2
Protein Name:	Cadherin EGF LAG seven-pass G-type receptor 2
Other Name:	CELSR2; CDHF10; EGFL2; KIAA0279; MEGF3; Cadherin EGF LAG seven-pass G-type receptor 2; Cadherin family member 10; Epidermal growth factor-like protein 2; EGF-like protein 2; Flamingo homolog 3; Multiple epidermal growth factor-like domains
Human Swiss Prot No.:	Q9HCU4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CDHF10 Polyclonal Antibody detects endogenous levels of CDHF10 protein.
Storage and Stability:	-20°C/1 year