

Antibodies Product Datasheet

Product Name:	CDHF11 Antibody
Catalog Number:	C12184
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CELSR3. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	CELSR3
Protein Name:	Cadherin EGF LAG seven-pass G-type receptor 3
Other Name:	CELSR3; CDHF11; EGFL1; FMI1; KIAA0812; MEGF2; Cadherin EGF LAG seven-pass G-type receptor 3; Cadherin family member 11; Epidermal growth factor-like protein 1; EGF-like protein 1; Flamingo homolog 1; hFmi1; Multiple epidermal growth factor-Q9NYQ7
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CDHF11 Polyclonal Antibody detects endogenous levels of CDHF11 protein.
Storage and Stability:	-20°C/1 year