

Antibodies Product Datasheet

Product Name:	c-FLIP Antibody
Catalog Number:	C12187
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CFLAR. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	CFLAR
Protein Name:	CASP8 and FADD-like apoptosis regulator
Other Name:	CFLAR; CASH; CASP8AP1; CLARP; MRIT; CASP8 and FADD-like apoptosis regulator; Caspase homolog; CASH; Caspase-eight-related protein; Casper; Caspase-like apoptosis regulatory protein; CLARP; Cellular FLICE-like inhibitory protein; c-FLIP; FAD
Human Swiss Prot No.:	O15519
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	c-FLIP Polyclonal Antibody detects endogenous levels of c-FLIP protein.
Storage and Stability:	-20°C/1 year