

Antibodies Product Datasheet

Product Name:	CIDE-3 Antibody
Catalog Number:	C12188
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CIDE. AA range:189-238
Concentration:	1 mg/ml
Gene Name:	CIDE
Protein Name:	Cell death activator CIDE-3
Other Name:	CIDE; FSP27; Cell death activator CIDE-3; Cell death-inducing DFFA-like effector protein C; Fat-specific protein FSP27 homolog
Human Swiss Prot No.:	Q96AQ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CIDE-3 Polyclonal Antibody detects endogenous levels of CIDE-3 protein.
Storage and Stability:	-20°C/1 year