

Antibodies Product Datasheet

Product Name:	CDK5 Activator-binding C48 Antibody
Catalog Number:	C12191
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDK5RAP2. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	CDK5RAP2
Protein Name:	CDK5 regulatory subunit-associated protein 2
Other Name:	CDK5RAP2; CEP215; KIAA1633; CDK5 regulatory subunit-associated protein 2; CDK5 activator-binding protein C48; Centrosome-associated protein 215
Human Swiss Prot No.:	Q96SN8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CDK5 Activator-binding C48 Polyclonal Antibody detects endogenous levels of CDK5 Activator-binding C48 protein.
Storage and Stability:	-20°C/1 year