

Antibodies Product Datasheet

Product Name:	HSF-27 Antibody
Catalog Number:	C12192
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDK5RAP3. AA range:341-390
Concentration:	1 mg/ml
Gene Name:	CDK5RAP3
Protein Name:	CDK5 regulatory subunit-associated protein 3
Other Name:	CDK5RAP3; IC53; MSTP016; OK/SW-cl.114; PP1553; CDK5 regulatory subunit-associated protein 3; CDK5 activator-binding protein C53; Protein HSF-27
Human Swiss Prot No.:	Q96JB5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSF-27 Polyclonal Antibody detects endogenous levels of HSF-27 protein.
Storage and Stability:	-20°C/1 year