

Antibodies Product Datasheet

Product Name:	Cofilin Antibody
Catalog Number:	C12219
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cofilin. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	CFL1/CFL2
Protein Name:	Cofilin-1/2
Other Name:	CFL1; CFL; Cofilin-1; 18 kDa phosphoprotein; p18; Cofilin; non-muscle isoform; CFL2; Cofilin-2; Cofilin, muscle isoform
Human Swiss Prot No.:	P23528/Q9Y281
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cofilin Polyclonal Antibody detects endogenous levels of Cofilin protein.
Storage and Stability:	-20°C/1 year