

Antibodies Product Datasheet

Product Name:	COX15 Antibody
Catalog Number:	C12230
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human COX15. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	COX15
Protein Name:	Cytochrome c oxidase assembly protein COX15 homolog
Other Name:	COX15; Cytochrome c oxidase assembly protein COX15 homolog
Human Swiss Prot No.:	Q7KZN9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	COX15 Polyclonal Antibody detects endogenous levels of COX15 protein.
Storage and Stability:	-20°C/1 year