

Antibodies Product Datasheet

Product Name:	CYP1A1/2 Antibody
Catalog Number:	C12249
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 1A1/2. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	CYP1A1/CYP1A2
Protein Name:	Cytochrome P450 1A1/2
Other Name:	CYP1A1; Cytochrome P450 1A1; CYPIA1; Cytochrome P450 form 6; Cytochrome P450-C; Cytochrome P450-P1; CYP1A2; Cytochrome P450 1A2; CYPIA2; Cytochrome P(3)450; Cytochrome P450 4; Cytochrome P450-P3
Human Swiss Prot No.:	P04798/P05177
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP1A1/2 Polyclonal Antibody detects endogenous levels of CYP1A1/2 protein.
Storage and Stability:	-20°C/1 year