

Antibodies Product Datasheet

Product Name:	CYP2A6 Antibody
Catalog Number:	C12258
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 2A6. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CYP2A6
Protein Name:	Cytochrome P450 2A6
Other Name:	CYP2A6; CYP2A3; Cytochrome P450 2A6; 1; 4-cineole 2-exo-monooxygenase; CYP1A6; Coumarin 7-hydroxylase; Cytochrome P450 IIA3; Cytochrome P450(I)
Human Swiss Prot No.:	P11509
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP2A6 Polyclonal Antibody detects endogenous levels of CYP2A6 protein.
Storage and Stability:	-20°C/1 year