

Antibodies Product Datasheet

Product Name:	CYP2C8 Antibody
Catalog Number:	C12262
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	IHC-p: 1/100 - 1/300. IF:1:50-200 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CYP2C8. AA range:241-290
Concentration:	1 mg/ml
Gene Name:	CYP2C8
Protein Name:	Cytochrome P450 2C8
Other Name:	CYP2C8; Cytochrome P450 2C8; CYP11C8; Cytochrome P450 IIC2; Cytochrome P450 MP-12; Cytochrome P450 MP-20; Cytochrome P450 form 1; S-mephenytoin 4-hydroxylase
Human Swiss Prot No.:	P10632
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP2C8 Polyclonal Antibody detects endogenous levels of CYP2C8 protein.
Storage and Stability:	-20°C/1 year