

Antibodies Product Datasheet

Product Name:	FGF-23 Antibody
Catalog Number:	C12369
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FGF23. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	FGF23
Protein Name:	Fibroblast growth factor 23
Other Name:	FGF23; HYPF; Fibroblast growth factor 23; FGF-23; Phosphatonin; Tumor-derived hypophosphatemia-inducing factor
Human Swiss Prot No.:	Q9GZV9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-23 Polyclonal Antibody detects endogenous levels of FGF-23 protein.
Storage and Stability:	-20°C/1 year