

Antibodies Product Datasheet

Product Name:	Gab 2/4 Antibody
Catalog Number:	C12388
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GAB2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	GAB2/GAB4
Protein Name:	GRB2-associated-binding protein 2/4
Other Name:	GAB4; GRB2-associated-binding protein 4; GRB2-associated binder 2-like; GAB2-like; GRB2-associated binder 4; GRB2-associated-binding protein 2-like; Growth factor receptor bound protein 2-associated protein 4; GAB2; KIAA0571; GRB2-associate
Human Swiss Prot No.:	Q2WGN9/Q9UQC2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gab 2/4 Polyclonal Antibody detects endogenous levels of Gab 2/4 protein.
Storage and Stability:	-20°C/1 year