

Antibodies Product Datasheet

Product Name:	GluR- δ 2 Antibody
Catalog Number:	C12395
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRID2. AA range:831-880
Concentration:	1 mg/ml
Gene Name:	GRID2
Protein Name:	Glutamate receptor delta-2 subunit
Other Name:	GRID2; GLURD2; Glutamate receptor delta-2 subunit; GluR delta-2 subunit
Human Swiss Prot No.:	O43424
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GluR- δ 2 Polyclonal Antibody detects endogenous levels of GluR- δ 2 protein.
Storage and Stability:	-20°C/1 year