

Antibodies Product Datasheet

Product Name:	eRF3a Antibody
Catalog Number:	C12400
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GSPT1. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	GSPT1
Protein Name:	Eukaryotic peptide chain release factor GTP-binding subunit ERF3A
Other Name:	GSPT1; ERF3A; Eukaryotic peptide chain release factor GTP-binding subunit ERF3A; Eukaryotic peptide chain release factor subunit 3a; eRF3a; G1 to S phase transition protein 1 homolog
Human Swiss Prot No.:	P15170
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eRF3a Polyclonal Antibody detects endogenous levels of eRF3a protein.
Storage and Stability:	-20°C/1 year