

Antibodies Product Datasheet

Product Name:	HSP70L Antibody
Catalog Number:	C12406
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HSPA1L. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	HSPA1L/HSPA1A
Protein Name:	Heat shock 70 kDa protein 1-like/Heat shock 70 kDa protein 1A/1B
Other Name:	HSPA1L; Heat shock 70 kDa protein 1-like; Heat shock 70 kDa protein 1L; Heat shock 70 kDa protein 1-Hom; HSP70-Hom; HSPA1A; HSPA1; HSPA1B; Heat shock 70 kDa protein 1A/1B; Heat shock 70 kDa protein 1/2; HSP70-1/HSP70-2; HSP70.1/HSP70.2
Human Swiss Prot No.:	P34931/P08107
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSP70 Polyclonal Antibody detects endogenous levels of HSP70 protein.
Storage and Stability:	-20°C/1 year