

Antibodies Product Datasheet

Product Name:	IL-20R β Antibody
Catalog Number:	C12439
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL20RB. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	IL20RB
Protein Name:	Interleukin-20 receptor subunit beta
Other Name:	IL20RB; DIRS1; Interleukin-20 receptor subunit beta; IL-20 receptor subunit beta; IL-20R-beta; IL-20RB; IL-20R2
Human Swiss Prot No.:	Q6UXL0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-20R β Polyclonal Antibody detects endogenous levels of IL-20R β protein.
Storage and Stability:	-20°C/1 year