

Antibodies Product Datasheet

Product Name:	IFI-56K Antibody
Catalog Number:	C12456
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IFIT1. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	IFIT1
Protein Name:	Interferon-induced protein with tetratricopeptide repeats 1
Other Name:	IFIT1; G10P1; IFI56; IFNAI1; ISG56; Interferon-induced protein with tetratricopeptide repeats 1; IFIT-1; Interferon-induced 56 kDa protein; IFI-56K; P56
Human Swiss Prot No.:	P09914
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IFI-56K Polyclonal Antibody detects endogenous levels of IFI-56K protein.
Storage and Stability:	-20°C/1 year