

## Antibodies Product Datasheet

|                               |                                                                                                                               |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | IFITM3 Antibody                                                                                                               |
| <b>Catalog Number:</b>        | C12462                                                                                                                        |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                             |
| <b>Applications:</b>          | WB, IF/ICC, ELISA                                                                                                             |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000.<br>ELISA: 1/40000. Not yet tested in other applications.    |
| <b>Modification:</b>          | Unmodified/Total                                                                                                              |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                             |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.         |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human IFM3. AA range:1-50                                 |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                       |
| <b>Gene Name:</b>             | IFITM3                                                                                                                        |
| <b>Protein Name:</b>          | Interferon-induced transmembrane protein 3                                                                                    |
| <b>Other Name:</b>            | IFITM3; Interferon-induced transmembrane protein 3; Dispanin subfamily A member 2b; DSPA2b; Interferon-inducible protein 1-8U |
| <b>Human Swiss Prot No.:</b>  | Q01628                                                                                                                        |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                       |
| <b>Specificity:</b>           | IFITM3 Polyclonal Antibody detects endogenous levels of IFITM3 protein.                                                       |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                  |