

Antibodies Product Datasheet

Product Name:	JNK2 Antibody
Catalog Number:	C12684
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAPK9. AA range:246-295
Concentration:	1 mg/ml
Gene Name:	MAPK9
Protein Name:	Mitogen-activated protein kinase 9
Other Name:	MAPK9; JNK2; PRKM9; SAPK1A; Mitogen-activated protein kinase 9; MAP kinase 9; MAPK 9; JNK-55; Stress-activated protein kinase 1a; SAPK1a; Stress-activated protein kinase JNK2; c-Jun N-terminal kinase 2
Human Swiss Prot No.:	P45984
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	JNK2 Polyclonal Antibody detects endogenous levels of JNK2 protein.
Storage and Stability:	-20°C/1 year