

Antibodies Product Datasheet

Product Name:	Cdc23 Antibody
Catalog Number:	C12818
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APC8. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	CDC23
Protein Name:	Cell division cycle protein 23 homolog
Other Name:	CDC23; ANAPC8; Cell division cycle protein 23 homolog; Anaphase-promoting complex subunit 8; APC8; Cyclosome subunit 8
Human Swiss Prot No.:	Q9UJX2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc23 Polyclonal Antibody detects endogenous levels of Cdc23 protein.
Storage and Stability:	-20°C/1 year