

Antibodies Product Datasheet

Product Name:	Beclin-1 Antibody
Catalog Number:	C13032
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Beclin-1 . at AA range: 110-190
Concentration:	1 mg/ml
Gene Name:	BECN1 GT197
Protein Name:	Beclin-1
Other Name:	beclin 1, autophagy related
Human Swiss Prot No.:	Q14457
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Beclin-1 Polyclonal Antibody detects endogenous levels of Beclin-1
Storage and Stability:	-20°C/1 year