

Antibodies Product Datasheet

Product Name:	Elongin A2/3 Antibody
Catalog Number:	C13048
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELOA2. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	TCEB3B/TCEB3C
Protein Name:	RNA polymerase II transcription factor SIII subunit A2/3
Other Name:	TCEB3B; TCEB3L; RNA polymerase II transcription factor SIII subunit A2; Elongin-A2; EloA2; Transcription elongation factor B polypeptide 3B; TCEB3C; TCEB3L2; RNA polymerase II transcription factor SIII subunit A3; Elongin-A3; EloA3; Transcr
Human Swiss Prot No.:	Q8IYF1/Q8NG57
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Elongin A2/3 Polyclonal Antibody detects endogenous levels of Elongin A2/3 protein.
Storage and Stability:	-20°C/1 year