

Antibodies Product Datasheet

Product Name:	FAF1 Antibody
Catalog Number:	C13054
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FAF1. AA range:531-580
Concentration:	1 mg/ml
Gene Name:	FAF1
Protein Name:	FAS-associated factor 1
Other Name:	FAF1; UBXD12; UBXN3A; CGI-03; FAS-associated factor 1; hFAF1; UBX domain-containing protein 12; UBX domain-containing protein 3A
Human Swiss Prot No.:	Q9UNN5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAF1 Polyclonal Antibody detects endogenous levels of FAF1 protein.
Storage and Stability:	-20°C/1 year