

Antibodies Product Datasheet

Product Name:	GADD 45γ Antibody
Catalog Number:	C13057
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GA45G. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	GADD45G
Protein Name:	Growth arrest and DNA damage-inducible protein GADD45 gamma
Other Name:	GADD45G; CR6; DDIT2; Growth arrest and DNA damage-inducible protein GADD45 gamma; Cytokine-responsive protein CR6; DNA damage-inducible transcript 2 protein; DDIT-2
Human Swiss Prot No.:	O95257
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GADD 45γ Polyclonal Antibody detects endogenous levels of GADD 45γ protein.
Storage and Stability:	-20°C/1 year