

Antibodies Product Datasheet

Product Name:	DAP-5 Antibody
Catalog Number:	C13062
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EIF4G2. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	EIF4G2
Protein Name:	Eukaryotic translation initiation factor 4 gamma 2
Other Name:	EIF4G2; DAP5; OK/SW-cl.75; Eukaryotic translation initiation factor 4 gamma 2; eIF-4-gamma 2; eIF-4G 2; eIF4G 2; Death-associated protein 5; DAP-5; p97
Human Swiss Prot No.:	P78344
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DAP-5 Polyclonal Antibody detects endogenous levels of DAP-5 protein.
Storage and Stability:	-20°C/1 year