

Antibodies Product Datasheet

Product Name:	Laminin α -1 Antibody
Catalog Number:	C13064
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAMA1. AA range:2501-2550
Concentration:	1 mg/ml
Gene Name:	LAMA1
Protein Name:	Laminin subunit alpha-1
Other Name:	LAMA1; LAMA; Laminin subunit alpha-1; Laminin A chain; Laminin-1 subunit alpha; Laminin-3 subunit alpha; S-laminin subunit alpha; S-LAM alpha
Human Swiss Prot No.:	P25391
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Laminin α -1 Polyclonal Antibody detects endogenous levels of Laminin α -1 protein.
Storage and Stability:	-20°C/1 year