

Antibodies Product Datasheet

Product Name:	Laminin α -3 Antibody
Catalog Number:	C13066
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAMA3. AA range:2571-2620
Concentration:	1 mg/ml
Gene Name:	LAMA3
Protein Name:	Laminin subunit alpha-3
Other Name:	LAMA3; LAMNA; Laminin subunit alpha-3; Epiligrin 170 kDa subunit; E170; Epiligrin subunit alpha; Kalinin subunit alpha; Laminin-5 subunit alpha; Laminin-6 subunit alpha; Laminin-7 subunit alpha; Nicein subunit alpha
Human Swiss Prot No.:	Q16787
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Laminin α -3 Polyclonal Antibody detects endogenous levels of Laminin α -3 protein.
Storage and Stability:	-20°C/1 year