

Antibodies Product Datasheet

Product Name:	Laminin α -5 Antibody
Catalog Number:	C13068
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAMA5. AA range:2381-2430
Concentration:	1 mg/ml
Gene Name:	LAMA5
Protein Name:	Laminin subunit alpha-5
Other Name:	LAMA5; KIAA0533; KIAA1907; Laminin subunit alpha-5; Laminin-10 subunit alpha; Laminin-11 subunit alpha; Laminin-15 subunit alpha
Human Swiss Prot No.:	O15230
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Laminin α -5 Polyclonal Antibody detects endogenous levels of Laminin α -5 protein.
Storage and Stability:	-20°C/1 year