

Antibodies Product Datasheet

Product Name:	Laminin β -2 Antibody
Catalog Number:	C13070
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAMB2. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	LAMB2
Protein Name:	Laminin subunit beta-2
Other Name:	LAMB2; LAMS; Laminin subunit beta-2; Laminin B1s chain; Laminin-11 subunit beta; Laminin-14 subunit beta; Laminin-15 subunit beta; Laminin-3 subunit beta; Laminin-4 subunit beta; Laminin-7 subunit beta; Laminin-9 subunit beta; S-laminin sub
Human Swiss Prot No.:	P55268
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Laminin β -2 Polyclonal Antibody detects endogenous levels of Laminin β -2 protein.
Storage and Stability:	-20°C/1 year