

Antibodies Product Datasheet

Product Name:	Laminin γ -3 Antibody
Catalog Number:	C13074
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAMC3. AA range:1361-1410
Concentration:	1 mg/ml
Gene Name:	LAMC3
Protein Name:	Laminin subunit gamma-3
Other Name:	LAMC3; Laminin subunit gamma-3; Laminin-12 subunit gamma; Laminin-14 subunit gamma; Laminin-15 subunit gamma
Human Swiss Prot No.:	Q9Y6N6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Laminin γ -3 Polyclonal Antibody detects endogenous levels of Laminin γ -3 protein.
Storage and Stability:	-20°C/1 year