

Antibodies Product Datasheet

Product Name:	Galectin-1 Antibody
Catalog Number:	C13075
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Galectin-1 . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	LGALS1
Protein Name:	Galectin-1
Other Name:	LGALS1; Galectin-1; Gal-1; 14 kDa laminin-binding protein; HLBP14; 14 kDa lectin; Beta-galactoside-binding lectin L-14-I; Galaptin; HBL; HPL; Lactose-binding lectin 1; Lectin galactoside-binding soluble 1; Putative MAPK-activating protein P
Human Swiss Prot No.:	P09382
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Galectin-1 Polyclonal Antibody detects endogenous levels of Galectin-1 protein.
Storage and Stability:	-20°C/1 year