

Antibodies Product Datasheet

Product Name:	Galectin-2 Antibody
Catalog Number:	C13076
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LEG2. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	LGALS2
Protein Name:	Galectin-2
Other Name:	LGALS2; Galectin-2; Gal-2; Beta-galactoside-binding lectin L-14-II; HL14; Lactose-binding lectin 2; S-Lac lectin 2
Human Swiss Prot No.:	P05162
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Galectin-2 Polyclonal Antibody detects endogenous levels of Galectin-2 protein.
Storage and Stability:	-20°C/1 year